

CES-520 Embedded Systems  
Fall 2006 – Lab #1

Assigned: August 22, 2006  
Due: August 29, 2006

*Lab: Getting Familiar with the Rabbit 3000*

1. Work through the RCM3000 Getting Started Manual.
2. Install the hardware and software.
3. Run the sample programs pertaining to the Rabbit 3000.
4. Find out what the RCM3000's MAC ID is.

Send an e-mail to the instructor that clearly states the MAC ID of your RCM3000 Development Kit. The MAC ID is six hex digits and is reported by one of the sample programs that you will run in this lab.